

# Fast System

Designed by Dr Joseph CHOUKROUN

&

**POLYSTERIBOX®**



**x1**  
Boite Fast System  
Box Fast System




**x2**  
Pince Kelly  
Kelly forceps



**x1**  
Pince Fast  
Fast Tweezer



**x1**  
Tournevis manuel  
Manual screwdriver



**x50**  
Etiquettes de stérilisation  
Sterilization labels



**x1**  
Plaque 12 trous  
Mesh 12 holes



**x2**  
Plaques 5 trous  
Meshes 5 holes



carré / squared L = 22 mm

**x1**  
Tournevis contre angle  
Screwdriver contra-angle



**x1**  
Foret  
Drill



3.5mm - 4mm - 5mm - 6mm -  
7mm - 8mm - 9mm - 11mm


**x19**  
Vis autotaraudantes  
Self-Tapping Screws  
3.5mm x2 - 4mm x2 - 5mm x2 - 6mm x2  
7mm x3 - 8mm x3 - 9mm x3 - 11mm x2



**x2**  
Coss  
Cross



**x1**  
Polysteribox S  
(208 x 111 x 54mm)



**x1**  
Tapis en silicone  
Silicone bed



**x100**  
Plots de sécurité  
Safety seals



**x1**  
Clé USB  
USB Stick

## La clé USB contient :

- Les instructions d'utilisation Fast system
- Les instructions d'utilisation Polysteribox
- Une vidéo de présentation Fast System
- Une vidéo de présentation Polysteribox

## The USB key contains:

- Fast system Instructions for use
- Polysteribox Instructions for use
- A Fast System video presentation
- A Polysteribox video presentation
- This product sheet

**SDI-TLB**  
TOUS LES BIOMATERIAUX.COM

49, rue Gioffredo - 06000 Nice, France  
Tel: 04 93 85 31 27 Mail : info@tlbm.fr